

# IG Kappa Gene Clonality Detection by PCR and Capillary Electrophoresis

## Clinical Indication and Relevance

- Aids in the diagnosis of B-cell malignancies
- Detection of disease or minimal residual disease monitoring

## Methodology

The assay is performed on isolated DNA with BIOMED-2 primers amplifying the VK, JK as well as intragenic and Kde region of the IG kappa gene. The gene rearrangements are detected by analyzing the PCR products using capillary gel electrophoresis.

## Sensitivity

The assay sensitivity for detection of clonal B-cell populations is 10% of lymphocytes.

## Turn-around Time

Five to seven working days

## Sample Requirements

### Collect

- Peripheral blood (PB): 3-5 mL, in purple top (sodium EDTA) tube; yellow top (ACD) tube acceptable.
- Bone marrow aspirate (BM): 1-3 mL, drawn into a syringe containing anticoagulant and then delivered in purple top tube (EDTA).
- Fresh or frozen tissue: fresh tissue should be obtained in a sterile manner, and a minimum 3 mm<sup>3</sup> of tissue is required. Put fresh tissues in culture medium or snap freeze (see transport section below).
- Formalin-fixed paraffin-embedded (FFPE) tissue blocks: send FFPE tissue blocks to the lab, or contact lab for instructions about cutting sections for molecular studies.

### Transport

- PB or BM samples should be delivered immediately to the lab at 2-8°C (wet ice or cold packs). PB and BM specimens should not be frozen.
- Fresh tissue samples should be delivered at room temperature in RPMI culture medium to the lab within 3 hours of collection, or snap frozen in liquid nitrogen at -70°C and packed in dry ice for delivery. Please do not allow frozen tissues to thaw.
- Formalin-fixed paraffin embedded (FFPE) tissue blocks can be delivered at room temperature.

### Stability

PB or BM samples: ambient - 1 hour; refrigerated - 48 hours.

### Unacceptable Samples

- Serum or plasma; frozen PB or BM; clotted blood; severely hemolyzed samples
- Tissue samples fixed in Zenker's, B5, or Bouin's fixatives
- Bone marrow biopsies decalcified in formic acid

### **CPT Code(s)**

- 81264: IGK@ (Immunoglobulin kappa light chain locus) (eg, leukemia and lymphoma, B-Cell) gene rearrangement analysis, evaluation to detect abnormal clonal population(s)
- G0452-26: Molecular pathology procedure; physician interpretation and report

### **References**

1. Evans PA et al. Leukemia. 21:207, 2007
2. Langerak AW et al. Leukemia. 26:2159, 2012
3. van Dongen JJ et al. Leukemia. 17:2257, 2003
4. van Krieken JH et al. Leukemia. 21:201, 2007